



中国认可  
国际互认  
检测  
TESTING  
CNAS L3457



苏州中纺联检验技术服务有限公司  
CNTAC Testing Service Co.,Ltd.(Suzhou)

地址：江苏省苏州市吴江区盛泽镇敦煌路 388 号汇赢金融大厦 1-2 楼  
Floor 1-2,Huiying Financial Building,No.388,Dunhuang Road,Shengze,Wujiang,Suzhou,Jiangsu.

<http://www.fcl-sz.org.cn>

Tel: +86 512 63081315 Fax: +86 512 63515926

## TESTING REPORT

Report security check web: [www.fcl-sz.org.cn](http://www.fcl-sz.org.cn)  
Security code:3431573636

Report No. ZFLJ0466319

Page 1 of 3



### Applicant Information (客户信息)

Applicant Name(委托单位) : Wujiang Tutaike Textile Finishing Co., Ltd.(吴江市涂泰克纺织后整理有限公司)  
Contact person(委托人) : Wu Haiyan(吴海燕)

### Sample Information (样品信息)

Sample Description(样品描述) : White protective clothing(白色隔离衣)  
Style Name(品名) : One-piece protective clothing(连体式隔离衣)  
Sample Quantity(样品数量) : One piece(一件)

- Sample receiving/Test started date (样品受理/检测起始日期) : 2020-03-27  
- Report date (报告日期) : 2020-03-27  
- The original sample is stucked on the last paper. (送检样品原样粘贴在本报告的末页)

### Test Performed (检验标准)

Judgement according to(评定依据):N/A  
- Selected test(s) as requested by applicants. (本实验室根据客户要求完成以下检测内容)

### Pronounce (声明)

The results shown in this report refer only to sample(s) tested unless otherwise stated.All the items are performed in the standard conditions, except the noted cases.(除非特别说明，测试结果只对来样负责。所有测试项目均在标准规定的环境下进行，有注明除外。)

Except for the requirement of the client, the test results and the conformity judgement of this report do not take the uncertainty of the test results into account.(除非客户要求，本报告检测结果及符合性判定不考虑测量结果的不确定度。)

Signed for and on behalf of

CNTAC Testing Center-Suzhou

Approved by (批准)	Main testing by (主检)	Prepare by (制单)
王楣艳	惠永久	钟冬民

### Pronounce(声明):

The results shown in this report refer only to sample(s) tested unless otherwise stated.(除非有其他说明，测试结果只对来样负责。)  
Test test report cannot be reproduced, except in full.(本报告只能全部复制，不得部分复制。)





中国认可  
国际互认  
检测  
TESTING  
CNAS L3457



苏州中纺联检验技术服务有限公司  
CNTAC Testing Service Co.,Ltd.(Suzhou)

地址：江苏省苏州市吴江区盛泽镇敦煌路 388 号汇赢金融大厦 1-2 楼  
Floor 1-2,Huiying Financial Building,No.388,Dunhuang Road,Shengze,Wujiang,Suzhou,Jiangsu.

<http://www.fcl-sz.org.cn>

Tel: +86 512 63081315 Fax: +86 512 63515926

<b>TESTING REPORT</b>	Report security check web: <a href="http://www.fcl-sz.org.cn">www.fcl-sz.org.cn</a> Security code:3431573636
Report No. ZFLJ0466319	Page 2 of 3



--- Test Result (测试结果) ---

1. Tensile Strength(拉伸强力) GB/T 3923.1-1997 Strip test(条样法)

5cm strip test, CRE tensile test machine, testing speed 100mm/min, gauge length 100mm.(5 厘米条样法, CRE 拉伸测试仪, 测试速度 100mm/min, 初始隔距 100mm)

Unit(单位): <N>

Lengthwise(长度方向)

Widthwise(宽度方向)

Breaking force(断裂强力)

137

81

Remark(备注):

The expired test standard for the test items used to be in the approved scope of CNAS& CMA, The test was conducted as client's request. (该项目所依据的检测方法为过期标准, 客户要求采用, 曾获 CNAS 认可及资质认定。)

Pronounce(声明):

The results shown in this report refer only to sample(s) tested unless otherwise stated.(除非有其他说明, 测试结果只对来样负责。)  
Test report cannot be reproduced, except in full.(本报告只能全部复制, 不得部分复制。)





中国认可  
国际互认  
检测  
TESTING  
CNAS L3457



苏州中纺联检验技术服务有限公司  
CNTAC Testing Service Co.,Ltd.(Suzhou)

地址：江苏省苏州市吴江区盛泽镇敦煌路 388 号汇赢金融大厦 1-2 楼  
Floor 1-2, Huiying Financial Building, No.388, Dunhuang Road, Shengze, Wujiang, Suzhou, Jiangsu.

<http://www.fcl-sz.org.cn>

Tel: +86 512 63081315 Fax: +86 512 63515926

**TESTING REPORT**

Report security check web: [www.fcl-sz.org.cn](http://www.fcl-sz.org.cn)

Security code:3431573636

Report No. ZFLJ0466319

Page 3 of 3



**Original (送检样品原样)**



=====**End of Report (检测报告结束)**=====

**Pronounce(声明):**

The results shown in this report refer only to sample(s) tested unless otherwise stated. (除非有其他说明, 测试结果只对来样负责。)  
Test report cannot be reproduced, except in full. (本报告只能全部复制, 不得部分复制。)

